

Blue Leather Cryogenic Gloves

Our blue leather gloves are available in both wrist and mid arm lengths. These two products are manufactured in a water-repellent leather, with ISO AIR felt and wool / cotton fleece lining.

These products have a highly efficient Porelle microporous membrane, which offers waterproof protection, but is also breathable.

Protection offered is down to -197° C.

Extra safety features have been added into the design of these Cryogenic blue leather gloves, as follows.

An additional safety feature is offered in the form of an adjustable strap on the cuff of the longer length gloves, to ensure complete protection against liquid nitrogen vapours. The full cuff length is 15 cm.

The wrist gloves have an added leather vein patch on the inner side of the knitted wrist cuff, giving added protection against cold environments.

Cryogenic Blue Leather Wrist Glove, Usage down to -197° C

| Size | Part Number | Price |
|----------------|-------------|-------|
| Medium, size 9 | 220245 | |
| Large, size 10 | 220246 | |

Cryogenic Blue Leather Glove, Usage down to -197° C

| Size | Part Number | Price |
|----------------------|-------------|-------|
| Medium, size 9 | 220248 | |
| Large, size 10 | 220249 | |
| Extra Large, size 11 | 220250 | |

1. UK Mainland carriage - 2 pairs of gloves

2. Quantity discount available.



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.