

Specification

Product name: H-ATDA - Adhesive crepe tapes

REF	Width x Length
H-ATDA	19 mm x 50 m

Materials and Properties

Backing:	Crepe paper
Indicator:	Steam indicator according to EN ISO 11140-1, class 1. Colour change from white to dark grey – black.
Colour:	beige
Adhesive:	normal adhesive
Total thickness:	0,165 mm
Adhesion to steel:	8 N / 25 mm
Tensile strength:	90 N / 25 mm
Elongation at break:	7 %
Toxicity:	toxic effects are not known

Compatibility:

Packaging material:	Crepe Paper	<input checked="" type="checkbox"/> suitable	<input type="checkbox"/> unsuitable
	Fleece	<input checked="" type="checkbox"/> suitable	<input type="checkbox"/> unsuitable
	Textiles material	<input type="checkbox"/> suitable	<input checked="" type="checkbox"/> unsuitable
Sterilization methods:	STEAM	<input checked="" type="checkbox"/> suitable	<input type="checkbox"/> unsuitable
	EO (Ethylene oxide)	<input type="checkbox"/> suitable	<input checked="" type="checkbox"/> unsuitable
	FORM (Formaldehyde)	<input type="checkbox"/> suitable	<input checked="" type="checkbox"/> unsuitable
	VH2O2 (Plasma)	<input type="checkbox"/> suitable	<input checked="" type="checkbox"/> unsuitable



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.