



**Dual-configuration for quick charge with
leak back absorbent protection
– avoid irritating hazardous shipment
supplemental charges!**

TW Dual Series Small Freezers

by IC Biomedical

- Universal canister compatibility
- Fits in standard shipping boxes and containers
- Rugged ribbed high-strength aluminum body, magneformed necktube design and more durable paint
- Superior vacuum performance with insulation technology that delivers maximum holding times
- Lockable lid with enlarged opening for seals or securing ties

Specifically designed for the every day requirements of laboratory and biorepository operations

SPECIFICATIONS	TW20 Dual	TW34 Dual	TW47-6 Dual	TW47-10 Dual
LN ₂ Capacity (L)	21	34	48	48
Static Hold Time (Days)	135	123	76	76
Vapor Hold Time (Days)	30	30	30	76
Evap. Rate (Liters/Day)	0.095	0.180	0.390	0.390
DIMENSIONS				
Necktube (in/mm)	2.2 / 55.9	3.5 / 88.9	5 / 127	5 / 127
Overall Dia. (in/mm)	14.5 / 368.3	18 / 457.2	20 / 508	20 / 508
Height (in/mm)	26 / 660.4	27 / 685.8	26.5 / 673.1	26.5 / 673.1
WEIGHT				
Empty (lb/kg)	26 / 11.8	35 / 15.9	41 / 18.6	41 / 18.6
Charged (lb/kg)	34 / 15.4	45 / 20.4	59 / 26.8	59 / 26.8
Full Weight (lb/kg)	63 / 28.6	94 / 42.6	120 / 54.4	120 / 54.4
CAPACITY				
Number of Canisters	6	6	6	10
Canister Diameter (in/mm)	1.65 / 41.9	2.8 / 71.2	4.0 / 101.6	2.8 / 71.1
.5 cc Straws on Canes	660	2100	3900	3200



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.