

## 1.0ml 96 round wells, round Bottom. Low profile



### PRODUCT SPECIFICATIONS

|                             |   |                      |                            |
|-----------------------------|---|----------------------|----------------------------|
| <b>Catalog Numbers:</b>     | IST-410-096MP<br>IST-410-096TP  | Pk of 50<br>Pk of 10 | NON-STERILE<br>NON-STERILE |
| <b>Product Description:</b> | <ul style="list-style-type: none"><li>- Working volume of 1.0ml</li><li>- Round well with 'Round bottom.</li><li>- High uniformity from well to well</li><li>- Raised well rims for reliable closing with heat sealing.</li><li>- Easy and reliable stacking</li><li>- Good centrifugation stability up to 6,000 × g for faster protocols and improved sample quality</li><li>- Manufactured under DNase/ RNase free environment without slip agents, plasticizers or biocides – Materials which could have a negative effect on bioassays</li><li>- Pure, virgin polypropylene guarantees good extractable performance, high resistance to chemicals, good mechanical stress and working with temperature extremes</li><li>- Autoclavable (121 °C, 20 min)</li></ul> |                      |                            |
| <b>Product Dimensions:</b>  | Length: 127.8mm<br>Width: 85.5mm<br>Height: 33mm  |                      |                            |
| <b>Packaging:</b>           | 10 plates per pack<br>10 packs per case (100 plates total)  |                      |                            |
| <b>Performance Testing:</b> | Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.  |                      |                            |



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.