



## Countertop Laboratory Essential Refrigerator 107L

GDPA code: **LSFSR107UK**

SKU code: **444410632**

- ✓ 2°C to 10°C temperature range
- ✓ Spark free interior
- ✓ External digital controller and temperature display
- ✓ Fitted lock and 2 keys
- ✓ Audible and visual temperature alarms



The Experts in  
Medical Refrigeration

## Specification: Countertop Laboratory Essential Refrigerator 107L

Gross Capacity (litres)	107
Net Capacity (litres)	103
Unit External Width (mm)	495
Unit External Depth (mm)	620
Unit External Height (mm)	830
Product Weight (kg)	31.5kg
Packed Width (mm)	560
Packed Depth (mm)	850
Packed Height (mm)	655
Packed Weight (kg)	35kg
Country of Origin	UK
Internal Dimensions (mm) (WxDxH)	403mm W x 700mm H x 400mm D
Temperature Range	2°C to 10°C
Colour	White/green
External Finish	Painted steel
Internal Finish	HIPS*
Type of plastic	SRL800 (HIPS/ABS mix)
Door Finish	Painted steel
Door Type	Solid
Door Style	Hinged reversible
Door Open Depth (mm)	1050mm
Replacement Door Seal	No
Open Door Sensor	Yes
Lock	Yes, 2 Keys
Adjustable Feet (number & position)	1 adjustable front foot
Porthole	No
Porthole Size	N/A
Porthole Location	N/A
Temperature Display on Panel	Yes
Temperature Alarm	Yes
Remote Alarm Contacts	Yes
Audible Alarms	Yes
Visual Alarms	Yes
Number of Shelves	2
Number of Shelf Positions	4
Shelf Size (mm) (WxDxH)	427mm W x 6mm H x 374mm D
Shelf Material	Wire
Number of Baskets	1
Basket Material	Wire
Basket Size (mm) (WxDxH)	—
Interior Light	No
Light Type	N/A
Light Activation	N/A
Fan	No
Spark Free	Yes
Auto or Manual Defrost	Auto
Controller Type	Digital, external
Number of Probes	1
Probe Type	Air
Climate Class	N
Refrigerant Type	R600a
CFC Free	Yes
HFC Free	No

**Specification:** Countertop Laboratory Essential Refrigerator 107L

<b>HCFC Free</b>	Yes
<b>Ammonia Free</b>	Yes
<b>UK 3 Pin Plug</b>	Yes
<b>Other Plug Types Available</b>	Yes
<b>Noise Level (dB)</b>	<43dB
<b>Energy Consumption (kwh)</b>	0.86 kW/24h
<b>Voltage Supply</b>	230
<b>Aux Output</b>	Yes
<b>Insulation Thickness (mm)</b>	35mm
<b>Insulation Type</b>	Cyclopentane
<b>Blowing Agent</b>	Yes
<b>Anti Microbial Handle</b>	No
<b>Optional Drawer System Available</b>	No
<b>Wall Brackets</b>	No
<b>Digital Lock</b>	No



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.