

SUFsg 5001

MediLine

Smart
Monitoring



Gross/net capacity	491/477 litre
Ext. dimensions in mm (w/d/h)	920/1082/1966
Int. dimensions in mm (w/d/h)	606/605/1300
Energy consumption in 365 days	2884 kWh
Ambient temperature range	+16°C to +32°C ¹
Refrigerant	R 290 / R 290 / R 170
Airborne acoustical noise emissions	47 dB(A)
Connection Rating	7.0 A
Voltage	230 V~
Cooling system	static
Defrosting	manual
Temperature range	-40°C to -86°C
Housing material/colour	steel/white
Door material	Steel
Liner Material	stainless steel
Type of Control	electronic controller
Temperature display	external digital
Fault: warning signal	visual and audible
Interface	LAN (Ethernet)
volt-free alarm contact	yes
SmartCoolingHub	prepared
Adjustable Shelving	3
Shelves, freezer compartment	4
Handle	Aluminium handle with integrated opening mechanism
Castors	castors with brake at the front, fixed castors at the rear
Lock	fitted
Door hinges	right
Insulation	103-180 mm
Defrost water drain	yes
Gross/net weight	308/254 kg
EAN-Nr.	9005382236550

¹ Indicates the ambient temperature at which the device can be operated safely.



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.