



Liquid Tight Sump Included

Use of caustic materials in the workplace can present many hazards, particularly if acids and alkalis are stored incorrectly or with other chemicals. These storage cabinets enable safe and identifiable storage in the workplace.

- “No snag” handles with 2-point locking.
- New regulation labels enclosed to be fitted by customer.
- Adjustable spill-retaining shelves for superior spill management*.
- Integral liquid-tight sump encourages the correct handling of spills through the use of absorbents*.

- Wall cupboards are pre-drilled for mounting (fixings not included).
- Stands raise cupboards by 543mm with an adjustable foot for uneven surfaces.
- Powder coated with white Germ Guard Active Technology Anti-Bacterial Paint.
- Reaction to Fire Classification EN 13501-1.
- It is not recommended to store acids and alkalis together in a single unit as corrosion may occur.
- Cupboard shelves have a UDL of 50kg.

*Except wall mounted cupboards.

Floor Cupboards - Dimensions/Codes

| Size H.W.D. | Shelves | Code | Price |
|---------------|---------|--------------|----------------|
| 1800.1200.460 | 3 | CA181246ZWXX | £571.34 |
| 1800.900.460 | 3 | CA189046ZWXX | £482.44 |
| 1200.900.460 | 2 | CA129046ZWXX | £417.24 |
| 900.900.460 | 1 | CA909046ZWXX | £318.86 |
| 900.460.460 | 1 | CA904646ZWXX | £248.92 |
| 700.900.460 | 1 | CA709046ZWXX | £292.78 |
| 700.350.300 | 1 | CA703530ZWXX | £208.62 |

Extra Shelves - Dimensions/Codes

| To Fit | Code | Price |
|----------------|--------------|---------------|
| W.1200 x D.460 | CAES1246ZWXX | £78.23 |
| W.900 x D.460 | CAES9046ZWXX | £67.56 |
| W.460 x D.460 | CAES4646ZWXX | £52.16 |
| W.350 x D.300 | CAES3530ZWXX | £41.49 |

Wall Cupboards - Dimensions/Codes

| Size H.W.D. | Shelves | Code | Price |
|-------------|---------|--------------|----------------|
| 570.850.255 | 1 | CA855725ZWXX | £277.37 |

Mobile Cupboards - Dimensions/Codes

| Size H.W.D. | Shelves | Code | Price |
|--------------|---------|--------------|----------------|
| 1040.900.460 | 1 | CA109046MWXX | £471.77 |
| 840.900.460 | 1 | CA849046MWXX | £445.69 |

Stands - Dimensions/Codes

| Size H.W.D. | Code | Price |
|-------------|--------------|----------------|
| 543.900.460 | CAST9046ZWXX | £119.72 |
| 543.460.460 | CAST4646ZWXX | £109.05 |
| 543.350.300 | CAST3530ZWXX | £100.75 |



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.