

## Swabbing Kits

High reach swabbing tools

A reproducible swabbing technique is vital for obtaining meaningful results during Cleaning Validation and environmental swabbing. Our swabbing kits allow the operator to make sure that a consistent area is swabbed.

**Swab Holder Kits** - include a Universal Swab Holder attached to a telescopic pole

**Swab Kit** - includes a Swabbing Template attached to a telescopic pole

### Swab Kit (3044B)



Each Swab Kit contains

- 1 x Telescopic Rod 3m
- 1 x Elbow connector
- 1 box Swabbing Templates\* (50 per box)

\* A choice of two template sizes is available. They are also available sterile or non-sterile.

The single use, plastic template will be in contact with the item being swabbed. In order to minimise the risk of sampling error the Swabbing Templates are made and packed within a Cleanroom.

- Made & Packed in a Cleanroom
- Individually Bagged
- Available sterile or non-sterile
- Full Batch Traceability
- 50 templates per box

### Swab Holder Kit (3044C)



Each Swab Holder Kit contains:

- 1 x Telescopic Rod 3m
- 1 x Elbow connector
- 1 x Universal Swab Holder (316L ss)

The Universal Swab Holder is designed to work with most makes of stick swabs and cloth swabs. The Swab Holder is made from hygienic 316L stainless steel.

Note: Swabs are not included with the Swab Holder Kit.



Swab Holder Kit  
(3044C-01)

Swab Kit 100  
(3044B-100)

Disposables



Swabbing  
Templates can  
be hand held

## Swabbing Kits

Part No.	Description
3044B-50	Swab Kit 50 - includes Swabbing Templates 50 x 50mm (50 per box)
3044B-50S	Swab Kit 50S - includes STERILE Swabbing Templates 50 x 50mm (50 per box)
3044B-100	Swab Kit 100 - includes Swabbing Templates 100 x 100mm (50 per box)
3044B-100S	Swab Kit 100S - includes STERILE Swabbing Templates 100 x 100mm (50 per box)
3044C-01	Swab Holder Kit

NOTE: Longer Telescopic Rods are available - ask for details

## SPARE SteriWare® Swabbing Templates

Part No.	Nominal Hole Dimensions	Material	Sterile	Box Quantity
8888H-50	50 x 50mm	HDPE	No	50
8888H-50S	50 x 50mm	HDPE	Yes	50
8888H-100	100 x 100mm	HDPE	No	50
8888H-100S	100 x 100mm	HDPE	Yes	50





# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.