

**1. Part Numbers:** 8126H-60A

**2. Photograph:**



### 3. Product Information:

Description	Wide Mouth HDPE Bottle 60ml Amber
Materials of Construction:	Bottle HDPE (conforms to FDA 177.1520 & USP Class VI) Cap PP (conforms to FDA 177.1520 & USP Class VI)
Manufacturing Environment:	Controlled environment within factory
Nominal Volume:	60ml
Brimful Volume:	68ml
Nominal Height:	85mm
Diameter:	39mm
Internal Neck Diameter:	21mm
Weight:	17g
Pack Size:	72 bottles
Method of Sterilisation:	N/A (these bottles are non-sterile)
BSE/TSE Free:	Yes
Supplied with lid screwed on:	Yes
Recommended Storage Conditions:	Dry and ambient temperature

### 4. Approval:

Signed:



Name: Kuldeep Flora  
 Position: Quality Manager  
 Date: 2-Jul-21



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.