

1. Part Number: 8180A-250

2. Photograph:



3. Product Information:

| | |
|---------------------------------|--|
| Description: | PP Tub 250ml |
| Use: | Wide mouth tub with screw on cap - for use with powders |
| Material of Construction: | Polypropylene (PP) – virgin <ul style="list-style-type: none"> • Conforms to FDA CFR 177.1520 • Conforms to EU Regulation 10/2011 • Conforms to EC Regulation 1935/2004 • BSE/TSE free |
| | White Colouring (1%) |
| Method of Construction: | Single piece, injection moulded |
| Nominal Volume: | 250ml |
| Brimful Volume: | 280ml |
| Method of Construction: | Single piece, injection moulded |
| Height: | 88mm |
| Lid Diameter: | 76mm |
| Nominal Weight: | 56g |
| Recommended Storage Conditions: | Dry and ambient temperature |

4. Approval:

Signed:



Name: Kuldeep Flora
 Position: Quality Manager
 Date: 6-Dec-21



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.