

**1. Part Numbers:** A125-5

**2. Photograph:**



**3. Product Information:**

Description	PTFE Bottle 5ml with screw on Cap
Material of Construction:	PTFE
Nominal Volume:	5ml
Overall Diameter:	22mm
Internal Diameter of Neck:	10mm
Overall Height of body:	32mm
Nominal Weight:	14g (body & lid)
Recommended Storage Conditions:	Dry and ambient temperature

**4. Approval:**

Signed:



Name: Kuldeep Flora  
Position: Quality Manager  
Date: 6-Dec-21



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.