


**1. Part Number:** A416-2

**2. Photograph:**



Dimensions

A = 119mm

B = 98mm

C = 228mm

D = 122mm

Wall Thickness  
= 1.2mm

### 3. Product Information:

Description:

Stainless Container 2L

Materials of Construction:

- Body: 316L stainless steel
- Lid: 316L stainless steel
- Clamp: 304 stainless steel
- Gasket: Silicone (FDA acceptable grade)

Method of Construction:

All welds ground and polished  
Crevice free interior

Surface Finish:

Better than 0.5 microns Ra

Nominal Volume

2 Litres

BSE/TSE statement:

All polishing compounds used in the manufacture of this product are of vegetable origin.

Weight:

2.0 Kg (body + lid + gasket + clamp)

Recommended Storage Conditions:

Dry and ambient temperature

Note: The neck of the container is fitted with a 4" diameter ferrule which conforms to BS4825-3

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 14-Aug-20



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.