


1. Part Number: A418-1

2. Photograph:




3. Product Information:

| | |
|-------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description | Stainless Container 1L with Internal GL45 Thread |
| Materials of Construction: | <ul style="list-style-type: none"> - Body: 316L stainless steel - Lid: 316L stainless steel - Clamp: 304 stainless steel - Gasket: Silicone (FDA acceptable grade) - Cap: 316L stainless steel - Cap Seal: Silicone (FDA acceptable grade) |
| Method of Construction: | All welds ground and polished Crevice free interior |
| Surface Finish: | Better than 0.5 microns Ra |
| Nominal Volume | 1 Litre |
| Dimensions | A = 119mm B = 101mm C = 152mm D = 104mm |
| Weight: | 1.6 Kg (body + lid + gasket + clamp) |
| Recommended Storage Conditions: | Dry and ambient temperature |
| Note: The neck of the container is fitted with a 4" diameter ferrule which conforms to BS4825-3 | |

4. Approval:

Signed:



Name: Kuldeep Flora
 Position: Quality Manager
 Date: 27-Oct-20



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.