

1. Part Number: 5005A-100

2. Photo:



3. Product Information:

Part No.	5005A-100
Description	Stainless Steel Dipper – 100ml
Material of Construction:	316L stainless steel
Method of Construction:	Fully Welded Construction
Surface Finish:	Cup better than 0.5 microns Ra Handle better than 1.0 microns Ra
Nominal Cup Volume:	100ml
Brimful Cup Volume:	100ml
Overall Length:	600mm
Max. Width (cup and handle)	49mm
Cup Diameter:	39mm
Cup Height:	100mm
Nominal Weight:	485g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin

4. Approval:

Signed:



Name: Kuldeep Flora
 Position: Quality Manager
 Date: 6-Nov-19



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.