

1. Part Number: A244-150**2. Photograph:****3. Product Information:**

Description:	Industrial Liquid Funnel
Material of Construction:	304 Stainless Steel
Top Diameter:	150mm
Outer Diameter of Spout:	20mm
Overall Length:	170mm
Nominal Spout Length:	85mm
Nominal Weight:	170g
Side Handle:	Spot welded to cone section of funnel
Recommended Storage Conditions:	Dry and ambient temperature

Note: Due the manufacturing method there is an internal crevice where the cone and spout are joined. These funnels are primarily designed for industrial applications.

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 30-Nov-21



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.