

1. Part Number: A304-450

2. Photograph:



3. Product Information:

Description	Economy Scoop 450ml
Nominal Volume:	450ml
Material of Construction:	304 stainless steel
Method of Construction:	Spot welded (NB. This method of construction will leave crevices)
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	255mm
Body Length:	145mm
Body Diameter:	75mm
Nominal Weight:	220g
Quantity per Pack:	1
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
 Position: Quality Manager
 Date: 8-Nov-19



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.