

	Product Specification	Specification Number	A504-2
		Revision	01

Part No. A504-2

Description Stainless Spoon

Nominal Volume: 10ml

Material of Construction: 304 stainless steel




Surface Finish: Better than 0.5 microns Ra

Overall Length: 210 mm

Bowl Length: 65 mm

Max. Bowl Width: 32 mm

Weight: 68 g

Specification Approved By:		Signature:	Name: Jonathan Fenner Position: Director Date: 22-May-14
			
Page 1 of 1			



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.