


1. Part Number: A736-1000

2. Photograph:



3. Product Information:

Description:	Stainless Beaker 1000ml
Method of Construction:	Crevice free, all welds ground and polished
Material of Construction:	316L Stainless Steel
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume:	1000ml
Dimensions:	
A - Internal Diameter:	99mm
B - Nominal Overall Diameter:	110mm
C - Nominal Height:	150mm
Wall Thickness	1.2mm
Nominal Weight:	480g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are free of animal products
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 30-Jun-21



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.