

	Product Specification	Specification Number	A764-30
		Revision	01

Part No. A764-30

Description Medicine Cup 30ml

Material of Construction: 304 stainless steel



Method of Construction: Single piece construction


Surface Finish: Better than 0.8 microns Ra

Top Diameter: 52mm

Base Diameter: 32mm

Height: 25mm

Nominal Weight: 14g

Specification Approved By:	Signature: 	Name: Kuldeep Flora Position: Quality Manager Date: 9-Jul-19
Page 1 of 1		



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.