

1. **Part Number:**
- 2000A-1000
  - 2000T-025
  - 2000T-050
  - 2000T-075
  - 2000T-100

2. **Photo:**



3. **Product Information:**

Description: U-D Cohesive Sampler & Tips  
Materials: 316L stainless steel  
(Note: the stainless steel spring is not grade 316L)  
Method of Construction: Crevice free, all welds ground and polished

**Sampler Body** (Part No: 2000A-1000)

Max. Diameter: 26mm  
Overall Length (sampler & tip): 1080mm  
External Finish: <1 Micron Ra  
Nominal Weight of Sampler (2000A-1000): 660g

**Tips**

Part No.	Length	Overall Diameter	Depth of Sample Chamber	Weight
2000T-025	30mm	26mm	3.2mm	60g
2000T-050	30mm	26mm	5.5mm	57g
2000T-075	30mm	26mm	7.6mm	56g
2000T-100	35mm	26mm	9.8mm	74g

BSE/TSE statement: All polishing compounds used in the manufacture of this sampler are of vegetable origin and so are BSE/TSE free.

Recommended Storage Conditions: Dry and ambient temperature

4. **Approval:**

Signed:



Name: Kuldeep Flora  
Position: Quality Manager  
Date: 7-Oct-20



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.