

# Sanitiser 9

## Ultrasonic Bath

Ultrawave's **Sanitiser-Series** offers outstanding performance for use in medical, healthcare, scientific, laboratory, chiropody, veterinary and many other professional environments. **Sanitiser-Series** Ultrasonic Cleaners are Class I Medical Devices, compliant to MDR (EU) 2017/745.

Cleans and disinfects, when used in conjunction with Ultraclean Disinfect+.

Effective against enveloped viruses (including Coronaviruses), bactericidal, mycobactericidal/tuberculocidal and fungicidal/yeastocidal.



The **Sanitiser 9** is a 9.5Litre heated, digital ultrasonic bath.

The menu driven control panel allows simple adjustment of cycle time (0 to 99 minutes), temperature (ambient to 70°C), ultrasonic power and frequency mode for precise and accurate processing.

Combining features such as degas mode and fluid level sensor, the **Sanitiser 9** offers repeatable cleaning performance with consistent cleaning results every time.



Specifications are subject to change for continual improvements

# Sanitiser 9

## Ultrasonic Bath

Construction Case	Powder coated 304 Stainless Steel
Construction Tank	304 Stainless Steel
Control	Digital LCD
Overall Dimensions	340 x 265 x 265mm (WxDxH)
Tank Dimensions	300 x 240 x 150mm (WxDxH)
Internal Basket Dimensions	265 x 210 x 110mm (WxDxH)
Maximum Capacity	10Litres
Working Capacity	9.5Litres
Ultrasonic Power	100W
Ultrasonic Power (per litre)	10.5W/Litre
Operating Frequency	32 - 38kHz <i>Frequency Leap</i>
Maximum Heating Capacity	Ambient - 70°C
Heater Power	150W
Timer	0 - 99 minutes
Mains Connection	230VAC
Drain Outlet	1/2" BSP valve (located at rear of machine)
Lid	ABS Plastic
Basket	304 Stainless Steel
Gross Weight (inc. lid, basket & packaging)	7.3kg
Net Weight	6.3kg
Warranty	1 Year as standard - extend to 3 Years if registered online with us (within 30days of purchase)



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.