

# LABORATORY BLENDERS

## MX1200XTX SERIES

### Features:

- Variable-speed control dial operates from 1,500 to 24,000 RPM
- PULSE feature overrides speed 10 at 30,000+ RPM
- Heavy-duty base with powerful, 3.75 peak input HP motor
- One-piece, dishwasher-safe, removable jar pad
- 2-liter, BPA-free, copolyester container
- Durable stainless steel drive coupling
- Also available with stainless steel containers

SKU	Technical Info	Approvals	Plug
MX1200XTX	120V, 60 Hz	cETLus, NSF	B ☺
MX1200XTXEE	230V, 50 Hz	CE, RoHS	F ☺
MX1200XTXEK	230V, 50 Hz	CE, RoHS	G ☺

*For additional models, see page 23.*





## ADDITIONAL MODELS

### One Liter Series – Single Speed

700BU	Base Only, 120V
800BU	Base Only, 230V (non CE)*
800EBU	Base Only, 230V (CE approved)*
800KBU	Base Only, 230V (CE approved)**

### One Liter Series - Two Speeds

7011BU	Base Only, 120V
7011HB	Base Only, Heavy Duty, 120V
7011HG	Heavy-Duty, Glass Container, 120V
7011HS	Heavy-Duty, Stainless Container, 120V
8011BU	Base Only, 230V (non CE)*
8011EB	Base Only, 230V (CE approved)*
8011KB	Base Only, 230V (CE approved)**

### One Liter Series, Two Speeds w/ Timer

7010BU	Base Only, 120V
7010HB	Base Only, Heavy Duty, 120V
7010HG	Heavy-Duty, Glass Container, 120V
7010HS	Heavy-Duty, Stainless Container, 120V
8010BU	Base Only, 230V (non CE)*
8010EB	Base Only, 230V (CE approved)*
8010KB	Base Only, 230V (CE approved)**

### CB15 Series

CB15BU	Base Only, 120V
CB15EB	Base Only, 230V*

### CB15V Series

CB15VXP	4-Liter Plastic Jar, 120V
CB15VXEB	Base Only, RPM Display, 230V*
CB15VXKB	Base Only, RPM Display, 230V**

### TBB145 Series

TBB145S6	2-Liter Stainless Container, 120V
TBB145S4	1.4-Liter Stainless Container, 120V
TBB145EB	Base Only, 230V*

### TBB160 Series

TBB160S6	2-Liter Stainless Container, 120V
TBB160S4	1.4-Liter Stainless Container, 120V
TBB160EB	Base Only, 230V*

### TBB175 Series

TBB175S6	2-Liter Stainless Container, 120V
TBB175S4	1.4-Liter Stainless Container, 120V
TBB175EB	Base Only, 230V*

### MX1000 Series

MX1000XTS	2-Liter Stainless Container, 120V
MX1000XTEES	2-Liter Stainless Container, 230V*
MX1000XTEKS	2-Liter Stainless Container, 230V**

### MX1100 Series



MX1100XTS	2-Liter Stainless Container, 120V
MX1100XTEES	2-Liter Stainless Container, 230V*
MX1100XTEKS	2-Liter Stainless Container, 230V**

### MX1200 Series

MX1200XTS	2-Liter Stainless Container, 120V
MX1200XTEES	2-Liter Stainless Container, 230V*
MX1200XTEKS	2-Liter Stainless Container, 230V**

### MX1300 Series

MX1300XTS	2-Liter Stainless Container, 230V
MX1300XTEES	2-Liter Stainless Container, 230V*
MX1300XTEKS	2-Liter Stainless Container, 230V**

\*F Plug Type  \*\*G Plug Type 



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.